
SIDE MOUNT BRACKETS FOR RECARO SPA/SPG/ProRacer HANS SEAT

The Brey-Krause R-9067 side mount bracket kit centers (1) Recaro SPA/SPG/Pro Racer HANS seat to the aftermarket Recaro sliders.

When mounting a seat in a 996 or Boxster it is necessary to use Brey-Krause seat track adapters to adapt Recaro to the chassis. Use part number R-9040 for Recaro part number 0.20.40 sliders.



Front of Car ----->

Use the supplied hardware to mount the side mounts to the seat track sliders.

- Mount the side rails to the seat track sliders.
- Mount the assembly to the seat.
- Place the seat, side mounts and track assembly into the car.
- Slide the seat back and mount the front holes of the sliders (finger tight).
- Slide the seat forward and mount the rear holes.
- Slide the seat back again and tighten the front bolts.

Performance driving events are by nature potentially hazardous to your health. Brey-Krause Mfg. Co. Inc. shall not be liable for injury, consequential, or other types of damages resulting from the use or misuse of our products. Proper installation of these parts is critical and should only be done by a professional shop.